

BOTA WALKER ALTA

WALKER AIR TALL

Our Walker Air Tall and Short offer an integrated pneumatic single bladder air system combined with circumferential compressive support and immobilization with a contoured-for-comfort liner.

Integrated single pneumatic bladder. Rocker bottom with ample tread for more natural ambulation. Extra heel width provides more surface contact and better stability. Padded insole for better shock absorption. Seamless foot liner for comfort. Breathable, moisture wicking, contoured-for-comfort liner. Steel reinforced rocker bottom.

Intended use:

Stable fracture of Talus, Tarsus, Tibia and Fibula, severe to moderate ankle sprain & strains, post-operative immobilization, metatarsal fractures.

Item No.	Product description	Size	Shoe Size
440421	Walker Air Tall	X-Small	Men's 2 - 4 / Women's 3 - 5 ½
440431	Walker Air Tall	Small	Men's 4 ½ - 7 / Women's 6 - 8
440441	Walker Air Tall	Medium	Men's 7 ½ - 10 ½ / Women's 8 ½ - 11 ½
440451	Walker Air Tall	Large	Men's 11 - 12 ½ / Women's 12 - 13 ½
440461	Walker Air Tall	X-Large	Men's 13+ / Women's 14+

