

Vission

Thumb Spica Laced

Anatomically contoured to “fit like a glove”

Our Vission Thumb Spica Laced provides a firm and comfortable support that is anatomically contoured to hold the wrist and thumb in a functional position to promote healing.

- Memory foam laminate for ultimate patient comfort
- Anatomically contoured for superior fit
- Malleable aluminum stays for customized fit
- Quick lace system
- Contoured, adjustable thumb-web strap
- Integrated cotton knit panel for ease of application



measurement: wrist circumference	size	catalog left	catalog right
5" - 5.75"	XS	300520L	300520R
5.75" - 6.25"	S	300530L	300530R
6.25" - 6.75"	M	300540L	300540R
6.75" - 7.75"	L	300550L	300550R
over 7.75"	XL	300560L	300560R