

Vission

Wrist Strapped 6" and 8"

Anatomically contoured to "fit like a glove"

Our Vission Wrist Strapped 6" or 8" provides firm and comfortable support. Both are anatomically contoured to hold the wrist in a functional position to promote healing.

- Memory foam laminate for ultimate patient comfort
- Anatomically contoured for superior fit
- Malleable aluminum stays for customized fit
- Easy to locate contrasting color D-rings
- Contoured, adjustable thumb web strap
- Integrated cotton knit panel for ease of application



measurement: wrist circumference	size	catalog left	catalog right
5" - 5.75"	XS	300020L	300020R
5.75" - 6.25"	S	300030L	300030R
6.25" - 6.75"	M	300040L	300040R
6.75" - 7.75"	L	300050L	300050R
over 7.75"	XL	300060L	300060R