



NeuroTherm[®]
SPINECATH[®]
Intradiscal Catheter



Trailblazing thermal technology

SPINECATH[♦] Intradiscal Catheter for Intradiscal Electrothermal Therapy (IDET[♦])



The catheter heating inside the disc



ELECTROTHERMAL[®] 20S Generator



Catheter inserted in disc under fluoroscopy

The SPINECATH Intradiscal Catheter delivers controlled levels of thermal energy during Intradiscal Electrothermal Therapy (IDET), an effective treatment for patients with chronic discogenic low back pain who have failed to respond to non-operative care. A clinically proven technique, Intradiscal Electrothermal Therapy (IDET), provides a unique minimally invasive method for delivering thermal energy to the disc annulus.

Minimally invasive spine therapy that hits the target—only from NeuroTherm. Our strategic intent is to be: The choice of surgeons worldwide for surgical techniques that reduce trauma and pain to the patient, reduce cost to the healthcare system, and provide better outcomes for surgeons. Please let us know how we can help.

Five-centimeter heating section for thermal conduction over a broad area

Temperature monitoring systems for controlled delivery of thermal energy

A 17-gauge introducer needle for unmatched access to the posterior annulus

Standard and extra-long lengths available to meet a variety of clinical needs

SPINECATH Intradiscal Catheter

REF	Description
7210440	SPINECATH Intradiscal Catheter
7210441	SPINECATH Intradiscal Catheter XL
7209601	Introducer Needle—Genesis II Standard 6" (box of 5)
7209600	Introducer Needle XL—Genesis II Standard 9"
7210644	ELECTROTHERMAL 20S Generator; includes: operations/service manual, 8-pin universal extension cable, power cord, and foot pedal

NeuroTherm[®]

30 Upton Drive, Suite 2
Wilmington, MA 01887, United States
1-888-655-3500
1-978-658-2378 FAX
usasales@neurotherm.com



429 Brighton Road
Croydon, Surrey CR2 6EU, United Kingdom
+44 (0) 20 8660 4374
+44 (0) 20 8660 9417 FAX
sales@neurotherm.com

Modemstraat 20
1033 RW Amsterdam, The Netherlands
+31-204-702-702
+31 204 702 701 FAX
www.nisales@neurotherm.com

www.neurotherm.com

NeuroTherm is a registered trademark of NeuroTherm, Inc.
© NeuroTherm Inc., March 2007

[♦]Trademark of NeuroTherm. Certain marks Reg. US Pat. & TM Off.



CE Mark indicates conformity with European Council Directive of 14 June 1993 concerning Medical Devices (93/42/EEC).

©2008 NeuroTherm, Inc.
07/08 1061268b